

Release Notes v9.1.5

We're constantly working to improve the experience in DWOS. Below you will find a summary of the items fixed or improved upon in this update.

Version compatible with the following Dental Wings products:

- 7Series lab scanners
- 3Series lab scanners
- iSeries lab scanners
- Standalone software licenses

BUG FIXES

- Fixed display issue of occlusion key scans on 7Series scanners
- Fixed display issue showing unrequired empty areas after removing scan faces
- Fixed screw hole display issue when merging custom abutments with waxing
- Fixed issue that could block the update of shared implant libraries
- Fixed issue affecting information confidentiality of cases sent to DWOS Connect
- Included the latest SDK update to support 3Dconnexion devices
- Improved margin detection for cases created with third-party scanners